




TAVDental
PRECISION. PASSION. PARTNERSHIP



CAD/CAM


User Manual



CAD/CAM in dentistry is used to improve the design and creation of dental restorations. Take advantage of TAV Dental CAD/CAM restoration products using Internal Hex platform or Conical connection.

Titanium Base Anti-Rotation, NP

Height (H)	NP
0.5mm	7200612
1.5mm	7200613
3mm	7200614
Abutment Screw, NP	7200358



Titanium Base Anti-Rotation, SP

Height (H)	SP
0.5mm	7200626
1.5mm	7200627
3mm	7200628
Abutment Screw, SP	7200086



Titanium Base Free-Rotation, NP

Height (H)	NP
0.5mm	7200615
1.5mm	7200616
3mm	7200617
Abutment Screw, NP	7200358



Titanium Base Free-Rotation, SP

Height (H)	SP
0.5mm	7200629
1.5mm	7200630
3mm	7200631
Abutment Screw, SP	7200086



Scan Body



Scan Body, NP 7200621
Abutment Screw, NP 7200358



Scan Body, SP 7200635
Abutment Screw, SP 7200086

Digital Analog



Model Implant Analog, NP
7200678



Model Implant Analog, SP
7200677

Titanium Base Anti-Rotation, SP

Height (H)	SP
1.5mm	7201301
3mm	7201302
Abutment Screw, SP	7201057



Titanium Base Anti-Rotation, WP

Height (H)	WP
1.5mm	7201314
3mm	7201315
Abutment Screw, WP	7201056



Titanium Base Free-Rotation, SP

Height (H)	SP
1.5mm	7201304
3mm	7201305
Abutment Screw, SP	7201057



Titanium Base Free-Rotation, WP

Height (H)	WP
1.5mm	7201317
3mm	7201318
Abutment Screw, WP	7201056



Scan Body



Scan Body, SP 7201309
Abutment Screw, SP 7201057



Scan Body, WP 7201322
Abutment Screw, WP 7201056

Digital Analog



Model Implant Analog, SP
7201334



Model Implant Analog, WP
7201335

Scan Body



Scan Body for Multi Unit, TD

7200636

Digital Analog



Model Multi Unit Analog, TD

7200682

Apply to Internal Hex and Connical Platform

CEMENTED		Ti Base for Multi-Unit, TD	7200433
		Multi-Unit Screw, TD	7200427
SCREW		Multi-Unit Bridge Screw, TD	7200434

*Use screw to connect final crown to the multi-unit abutment.

Scan Bodies

TAV Dental Scan bodies are designed to indicate the position of dental implant (chairside) or analog (lab) in CAD/CAM scanning procedures.

Advantages

- ▶ Unique surface finish combined with special coating for superior scanning
- ▶ Reusable - Cost effective
- ▶ Laser marks for easy platform identification
- ▶ Short surface for indication with implant hex
- ▶ Unique non symmetric geometry for easy scan



SP - Standard Platform
IH - Internal Hex.

step 1

Before placing the scan body ensure that all components are clean and sterilized.

step 2

Make sure to select the right scan body according to the implant connection and platform.

step 3

Place the scan body in the dental implant or Analog and verify proper fit.

step 4

Insert the fixation screw using screw driver key and hand-tighten to maximum 15N·cm

step 5

If a single-tooth restoration is planned, ensure that the angled surface of the scan body is oriented buccally.

step 6

Perform the scan procedure according to the scanner in use and send files to lab (once performed chairside).

step 7

Following the scanning, remove the scan body and perform the sterilization according to IFU.

TAV Dental CAD/CAM libraries includes leading CAD/CAM providers.

Download and extract library files enter TAVDental link:
www.tavdental.com

The logo for exocad, featuring the word "exocad" in a bold, dark blue, lowercase sans-serif font. The logo is enclosed in a thin blue rectangular border.The logo for 3shape, featuring the word "3shape" in a dark blue, lowercase sans-serif font, followed by a red triangle icon. The logo is enclosed in a thin blue rectangular border.

Titanium Base

- ▶ Anodized Titanium
- ▶ Suitable for cement or screw retained restoration
- ▶ For single tooth (anti rotation) and bridge (free rotation) restoration
- ▶ Multi-Unit restoration
- ▶ Wide bonding surface for reliable adhesion
- ▶ Apply to internal hex and conical platform



Multi-Unit

Available Platform | all*



Ti-Base Internal Hex.

Available Platform | NP, SP



Ti-Base Conical Connection

Available Platform | SP, WP



* Apply to internal hex and conical platform

www.tavdental.com