



When Quality meets Prestige

TAV Dental Presents: iCone.

The new conical connection implant system. An advanced combination of functionality and uncompromising quality.

Strong and stable fit between implant and Abutment ♦ 11° cone that seals off the interior of the implant, minimizing micromovements and microleakage ♦ Internal Double Hexagon – features 12 indexing positions ♦ Platform switching design ♦ Available in a versatile range of diameters and lengths for all indications ♦ Color coded platforms for easy identification ♦ Restorative components for long-term function and esthetics.



www.TAVDental.com
info@tavdental.com
Tel. 04-9808615



TAVDental
PRECISION. PASSION. PARTNERSHIP